

# THE COMMUNICATOR

Introducing the world's first FDA-approved medical mask with a clear window. ASTM F2100 Level 1 Surgical Facemask.

## Wearing The Communicator

1. The plastic goes to the outside; the non-woven goes to the inside to absorb moisture from breathing.
2. The nose clip goes to the top, on the nose, which ensures a snug fit.
3. Be sure to pull the lower portion of the mask over the chin so the window will sit where the mouth area is visible, thereby enabling others to see facial expressions more easily.



[www.safenclear.com](http://www.safenclear.com)